

SANITARY DIVISION OF  
**YANGA**<sup>®</sup>  
ITALIAN TEXTILE GROUP



MANUFACTURED BY YANGA  
ITALY



COMPANY  
WITH QUALITY SYSTEM  
CERTIFIED BY DNV  
= ISO 9001/2008 =

**YANGA SRL**



VIA MONTE ORFANO 21

13864 CREVACUORE [BIELLA] ITALY



WEBSITE [www.yanga.it](http://www.yanga.it)




TEL +39 015 768256



FAX +39 015 768910



E-MAIL [yanga@yanga.it](mailto:yanga@yanga.it)



**REDIHIP**<sup>®</sup>

**YANGA**<sup>®</sup>

BENDAGGIO POST-OPERATORIO PER L'ANCA

HIP POST-OPERATIVE BANDAGE

Bendaggio post-operatorio che assicura sostegno e protegge le medicazioni dopo interventi all'anca, con attacchi rapidi in Velcro®. Prodotta con speciale tessuto morbido di cotone elasticizzato che garantisce ottima traspirabilità e aderisce perfettamente al corpo.

Può anche essere usata per mantenere "in loco" le medicazioni (ad es. dopo prelievo di midollo osseo).

La chiusura in velcro consente al paziente di applicare correttamente la cintura anche senza assistenza.

**HA UNA TRIPLICE AZIONE:** Contenitiva, sostenitiva, termica.

Disponibile in 4 taglie. Reversibile per uso su entrambe le anche.

#### VANTAGGI E BENEFICI:

- ⊙ Applicazione rapida e sicura.
- ⊙ Mantiene la medicazione in sede.
- ⊙ Aiuta il drenaggio delle ferite.
- ⊙ Riduce il rischio di ematoma.
- ⊙ Permette i movimenti delle articolazioni.
- ⊙ Garantisce comodità e sicurezza al paziente.
- ⊙ Consente un agevole accesso ai servizi igienici.
- ⊙ Lavabile e autoclavabile.



SIZE	HIP CIRC. / CIRC. ANCA	THIGH CIRC. / CIRC. COSCIA	CODE / CODICE PRODOTTO
S	75 - 90 cm	60 cm	HIP-S
M	90 - 110 cm	70 cm	HIP-M
L	110 - 135 cm	80 cm	HIP-L
XL	135 - 165 cm	90 cm	HIP-XL

Post-operative bandage that provides support and protects the hip region following injuries of surgery. with Velcro® fastening for easy fixation. Made of special soft stretch knitted cotton that ensures excellent breathability and adheres perfectly to the body. It can also be used to maintain in situ dressing (for example, after the offspring of bone marrow). Velcro® fastenings allow the patient to correctly apply the belt without assistance.

**IT HAS A TRIPLE ACTION:** containment, contenitive and thermal.

Available in 4 sizes and reversible for use on either hip.

#### ADVANTAGES AND BENEFITS:

- ⊙ Quick secure application.
- ⊙ Holds dressing in place.
- ⊙ Assists wound drainage.
- ⊙ Reduces the risk of haematoma.
- ⊙ Allows joint movement.
- ⊙ Provides comfort and security for the patient.
- ⊙ Allows easy access to toilet area.
- ⊙ Washable and autoclavable.

